EN DE
FR ES

MultiLab[®] Pro IDS 20

MultiLab[®] Pro IDS 40

Digital meter for IDS sensors

en	Quick Start Guide.....	3	fr	Guide de démarrage rapide.....	11
de	Schnellstart-Anleitung.....	6	es	Guía de inicio rápido.....	15

1 About this quick start guide

The purpose of this quick start guide is to give basic information on the installation and operation of the device. For more information, see the installation, operation, and maintenance manual. The installation, operation, and maintenance manual is accessible over the QR code that is given on this document.

2 Safety

2.1 Safety terminology

About safety messages

Read and follow the safety messages and rules before handling the product. They are published to help prevent these hazards:

- Personal accidents and health problems
- Damage to the product and its surroundings
- Product malfunction

2.2 User safety

Introduction

All government regulations, local health, and safety directives must be followed.

Prevent danger due to electricity

WARNING: Electrical Hazard

There is a risk of electrical shock or explosion if the electrical connections are not correctly carried out, or if there is fault or damage on the product. Visually inspect equipment for damaged cables, cracked casings or other signs of damage. Make sure that electrical connections have been correctly made.

All danger due to electricity must be avoided. Electrical connections must always be carried out in compliance with the following:

- The standard connections shown in the product documentation that is delivered together with the product.
- All international, national, state, and local regulations. For details, consult the regulations of local electricity suppliers.

2.3 Authorized use

The usage of the device must satisfy the following conditions:

- Always measure the pH, oxidation-reduction potential (ORP), conductivity, dissolved oxygen

(DO), and turbidity only in a laboratory environment.

- The device operation must adhere to the instructions and technical specifications that are given in this manual.
- Always obey the safety instructions of the sensors. For more information, see <https://www.xylenanalytics.com>.

2.4 Unauthorized use

Do not operate the device in the following cases:

- A visible damage is present on the device or the power adapter.
- The device was stored under unspecified conditions for an extended period.

Table 1: Storage temperature

Storage temperature	Value
Temperature range	-25 °C to +65 °C (-13 °F to +149 °F)

2.5 Data security

To protect the digital information against an unauthorized access, damage, or theft, obey the following points.

- The device must operate in a segmented network environment that is isolated from external internet access.
- The access must be strictly limited to authorized personnel or approved teams that are enforced through network controls and authentication systems.

Restrictions:

- No inbound or outbound internet connectivity is allowed.
- The device must not be exposed to any publicly routable IP address.
- All inter-system communications must be tightly controlled, auditable, and limited to predefined ports and protocols.

2.6 Requirements for safe operation

For the safe operation adhere to the following conditions:

1. Use the device only for the approved measurements.
2. Operate the device under the given environment conditions given in the following table.

Item	Value
Operating temperature	0 °C to +40 °C (32 °F to +104 °F)
Operating humidity	Average per year ≤ 80%

3. Never bend the adapter cable.

Bending can damage the electrical wires which leads to a fire formation.

2.7 Warranty

The warranty expires in the following cases:

- Use of a power adapter, other than the recommended one.
- Any physical damage on the device or the power adapter.
- Try to open the device.


3 Product description

3.1 Device layout

EN **Figure 1: MultiLab Pro IDS 40 touchscreen**

Figure 2: MultiLab Pro IDS 40 socket field

1. IDS sensor ports, A-D
2. Network port
3. USB-A port, 2x
4. Power port USB-C
5. Power switch

4 Set up the device

4.1 Examine the delivery

Examine the package

1. Examine the package for damaged or missing items upon delivery.
2. Record any damaged or missing items on the receipt and freight bill.
3. If anything is out of order, then file a claim with the shipping company.

If the product has been picked up at a distributor, make a claim directly to the distributor.

Examine the unit

1. Remove packing materials from the product.
Dispose of all packing materials in accordance with local regulations.
2. To determine whether any parts have been damaged or are missing, examine the product.
3. If there is any issue, then contact a sales partner.

4.2 Turn on the device

NOTICE:

- Use only the original power adapter from the package or a compatible power bank with at least 15 W
- Check that the power adapter is in good condition.

Parameter	Value
AC input	100 V-240 V, 50 Hz-60 Hz
DC output	5 V, 3 A, 15 W

1. Connect the correct primary plug to the power adapter.
There are four primary plugs that are supplied in the package.
2. Connect the power adapter cable to the power port USB-C on the device.
3. Connect the power adapter to the power outlet.
4. On the device, press the power switch.
The home screen is shown.

4.3 Principles of operation

- All Xylem blue elements are active elements, for example the edit icon . The user can execute the function or edit the setting.
- If the device is operated with an active user administration, then some elements are not accessible based on the permissions for the user.
- Active elements in the user interface can be menus, buttons, and icons.
 - Icons are shown as images, for example, .
 - Texts in the user interface, for example buttons or menu names, are displayed as **bold text** in the manual.

Here is an example of how the navigation is presented in this manual:

Short navigation:

Go to  > **System** > **Device settings** > **Language**.

Long navigation:

1. Tap the icon  in the status bar.
2. Tap the **System** button on the home screen.

- Tap the **Device settings** button on the system screen.
- Tap the **Language** button in the menu bar.

Touch gestures and actions

Table 2: Touch gestures and actions

Gesture	Description
Tap	The tap gesture selects an item or option on the display.
Swipe	The swipe gesture from the top, the left or the right opens the mostly used screens or menu options.

Elements on the home screen

Figure 3: Elements on the home screen


Table 3: Elements on the home screen

Position	Element	Description
1	Measurement preview	The preview of the measured values from all connected sensors is shown.
2	Button	The home screen has predefined buttons.
3	Top swipe indicator	The swipe indicator shows that a swipe gesture from the top edge is supported. If the swipe gesture is executed, then the settings screen is shown.
4	Right swipe indicator	The swipe indicator shows that a swipe gesture from the right edge is supported. If the swipe gesture is executed, then the database screen is shown.
5	Status bar	For more information, see Elements in the status bar on page 5.
6	▶	This button opens the measurement screen. The selected measurement method is active.
7	Left swipe indicator	The swipe indicator shows that a swipe gesture from the left edge is supported. If the swipe gesture is executed, then the Calibration methods menu is shown.

Elements in the status bar

The status bar displays icons for frequently used functions and status information for a quick overview. Functions are displayed in Xylem blue and orange.

Figure 4: Elements in the status bar



Table 4: Elements in the status bar

Position	Element	Description
1	Home icon	The icon opens the home screen.
2	Log out icon	The icon logs out the user without confirmation. The icon is gray when the function User Management is not active.
3	User icon	The icon opens the user settings to change the password and see the permissions. The icon is gray, when user management is not active.
4	User name	The username of the logged-in user is shown.
5	Signal strength icon	The icon opens the menu for the management of wireless sensors. For more information, see Connect a wireless sensor . If wireless sensors are connected, then the number of sensors is shown next to the icon.
6	Battery icon	The icon opens the installation, operation, and maintenance manual.
7	Date and time	The current date and time is shown.

4.4 Change the device settings

- Go to > **System** > **Device settings** > **Language**. The list of languages is shown.

Figure 5: Device settings

- Select the language.

3. Tap **OK**.

The language is changed.
The system screen is shown.

4. Go to **Device settings > Date / Time**.

5. Turn on the toggle switch to set the date and time manually.

6. Select the date and time.

7. Tap **OK**.

The date and time are set.

4.5 Connect the IDS sensors

The device supports multiple IDS sensors:

Model	Number of sensors	Sensor connection type
MultiLab Pro IDS 20	Up to 2	• Wired • Wireless
MultiLab Pro IDS 40	Up to 4	

1. Connect wired sensors directly to the IDS sensor ports A-D on the device.

2. Connect wireless sensors with the IDS WLM system accessory.

For more information, see the IDS WLM system manual.


The IDS sensor is shown in the measurement preview.

4.6 Access the measurement screen

On the home screen, tap the  button.

The measurement screen is shown.

Figure 6: Measurement screen - MultiLab Pro IDS 40


4.7 Start a calibration

1. Go to  > **CAL**.

The **Calibration methods** menu is shown.

Figure 7: Calibration methods menu

1. **Calibration** button
2. Start icon
3. Edit icon

2. Tap the  icon next to the name of the calibration method to be used.

The **Calibration settings** menu for the selected sensor is shown.

3. If necessary, change the settings.

4. Apply changes or cancel.

5. Tap the  icon next to the name of the calibration method to be used.

– If there is only one sensor connected that can be calibrated with that method, then the summary of the calibration method is shown.

– If there are more sensors connected that can be calibrated with that method, then the list of these sensors is shown. Select the name of the sensor to be calibrated.

A summary of the calibration method is shown.

6. Tap **Start**.

The procedure of the selected calibration method starts. Operating instructions for running the method are shown on the display.

7. Follow the instructions that are shown on the screen to complete the calibration.

For more information, see the installation, operation, and maintenance manual.

4.8 Show the installation, operation, and maintenance manual

Tap the  icon in the status bar.

The installation, operation, and maintenance manual are shown.

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© 2025 Xylem Inc. WTW is a trademark of Xylem Inc. or one of its subsidiaries. All
other trademarks or registered trademarks are property of their respective owners.
Xylem Analytics Germany GmbH +49 881 183-325